


## Catalog Number

9509S-10C-B

**Caps, Screw Top**, with Silicone Rubber / PTFE, ultra-pure Septa, fitted in "knurled", polypropylene, Light Blue Caps. For use with 9-425mm thread, 12x32mm autosampler vials. AQR Brand. 100/PK.



### Additional Info

When highest purity in a Septa is required. No color added to the rubber or Teflon for improved purity. Good choice for LCMS, GCMS & LCMS/MS.

*For Waters® Alliance, Acquity™ and all other autosamplers that require a Pre-Slit Septa use catalog# [9509S-30C-B](#).*

Click [HERE](#) for case quantities.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$0.00**

Printed On 2024/05/20