

## Catalog Number

9509S-10C-KSR

**Caps, Screw Top**, Solid, with an ultra-pure, Silicone Rubber / PTFE liner, in a closed Black Cap. For use with 9-425mm thread, 12x32mm autosampler vials. AQR Brand 100/PK.



### Additional Info

Ultra-pure liner, No color dyes added to the Rubber or PTFE. Good choice for storage when inertness and a good seal is required. Silicone rubber is not colored.

Turns 12x32mm, 2ml autosampler vials into storage bottles.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$0.00**

Printed On 2024/05/20