

Catalog Number

9509S-30C-B-M

Caps, Screw Top, with Silicone Rubber / PTFE, ultra-pure, Pre-Slit Septa, fitted in "knurled", polypropylene, Light Blue Caps. For use with 9-425mm thread, 12x32mm autosampler vials. AQR Brand. 1000/CS.



Additional Info

Highest purity Septa on the market, with no color added to the Rubber or Teflon for improved purity. Good choice for LCMS & LCMS/MS.

Pre-slit for easy penetration by autosampler needles and recommended for Waters® Alliance, Acquity™ and all other autosamplers that require Pre-Slit Septa. For lower cost, non-slit version please click [HERE](#) for more information.

Click [HERE](#) for smaller pack quantities.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$0.00

Printed On 2024/05/20