

## Catalog Number

9509S-WAV-RS

**Vials, Screw Top, Glass. Amber, 2ml, with Write-On Patch and fill lines, 9-425mm threads, 12x32mm outer dimensions. For use as an autosampler vial. RSA Brand. 100/PK.**



### Additional Info

First hydrolytic, Class A, Type N51A, Borosilicate Glass, amberized. Write-On Patch and fill lines are ceramic and will resist smearing or scratching when working with solvents.

Bottom depths and heights are designed, and QC controlled for maximum recovery and needle safety. Low adsorbing and will not add Na to sample solutions or compounds that cause pH changes.

Click [HERE](#) for more information about RSA autosampler vials.

*This is our most popular 2ml amber autosampler vial because of the consistency of the glass chemistry from batch to batch, vial to vial.. Reduced surface activity is tested for each lot and certified. Excellent for low abundance analytes, pH sensitive analytes and for LCMS.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$0.00**

Printed On 2024/05/20