

## Catalog Number

9509S-WCV-AQ

**Vials, Screw Top, Glass.** Clear, 2ml, with Write-On Patch and fill lines, 9-425mm threads, 12x32mm outer dimensions. For use as an autosampler vial or storage bottle. AQ Brand. 100/PK.



### Additional Info

First hydrolytic, Class A, Type-33 Borosilicate Glass. Write-On Patch and fill lines are ceramic and will resist smearing or scratching when working with solvents.

Compares to HPLC "Certified" vials for performance but lower priced and a great value. Bottom depths and heights are designed, and QC controlled for maximum recovery and needle safety.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$37.85**

Printed On 2024/05/20