

Catalog Number

9512C-0-0DCV

Vials, Snap Top, Glass. Clear, 1.8ml, MRQ Style. Deactivated. An 11mm snap-ring and 12x32mm outer dimensions. For use as an autosampler vial. MicroSolv Brand. 100/PK.



Additional Info

First hydrolytic, Class A, Type-N51A, Deactivated / Silanized, Borosilicate Glass. The MRQ Style Vials have a fill volume of 1.2ml and a very low residual volume of 5ul due to the center-draining effects and the conical bottom with a well that is slightly larger than autosampler vial needles. Side entry port needles will recover the maximum amount during use.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$267.85

Printed On 2024/05/20