

Catalog Number

9512C-0-0DCV-A

Vials, Crimp Top, Glass. Clear, 1.8ml, Max Recovery™. Deactivated. An 11mm Serum Finish with 12x32mm outer dimensions. For use as an autosampler vial. AQ Brand. 100/PK.



Additional Info

First hydrolytic, Class A, Type 33, Deactivated / Silanized Borosilicate Glass. Features an 8ul residual volume for most autosamplers due to the center draining interior, smooth walls (no ledges) and the conical bottom. This product is manufactured in an ISO 9001:2015 Certified facility.

Compares to Agilent Maximum Recovery vials.

Use with 11mm aluminum crimp caps. These lower cost center draining vials have serum finish tops that provide the best seal for chromatography, when used properly with correct aluminum crimp caps.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$346.98

Printed On 2024/05/20