

Catalog Number

9512S-0CV-P

Vials, Screw Top, Glass. Clear, 1.5ml, U Shape. 9-425mm threads, 12x32mm outer dimensions and a fluted stand. For use as an autosampler vial. MicroSolv Brand. 100/PK.



Additional Info

First hydrolytic, Class A, Type N51A, Borosilicate Glass.

Compares to Champagne Style Vials. Features a 15ul residual volume for most autosamplers due to the center draining, inner U shape.

We recommend to use these vials with [Soft-Guard™](#) 9-425mm Screw Caps.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$86.84

Printed On 2024/05/20