

## Catalog Number

95296-05N

**Inserts, Low Volume**, Glass. Clear, 350ul with a wide, conical point and 6x29mm outer dimensions for use in the Micro-Sample Management System. Replacement Inserts. U-2D Brand. 100/PK.



### Additional Info

Inserts are 1st Hydrolytic, Class A, TYPE 33 Borosilicate Glass made for use with the U-2D Micro-Sample Management System Rack or with 12x32mm, 2ml AQ Brand autosampler vials.

Wide conical points help with fast filling and minimizes bubble formation in the insert.

Also available pre-loaded in U-2D Brand racks of 96; click [HERE](#) for more information.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$125.11**

Printed On 2024/05/20