

## Catalog Number

95296-1mL-3R

**Micro-Sample Management System.** This set includes 96, 1000ul Precision, Glass Inserts with wide conical points in a Rack with a Base and a Pre-Slit, Solvent Resistant, Polymeric Plug Style Sealing Mat. U-



### Additional Info

This Set is a convenient way to purchase all that is needed to manage your Micro-Samples or use this set in a 96 Well Plate configuration. 1ml Inserts are 1st Hydrolytic, Class A Type 33 Borosilicate Glass. Use with most HPLC compatible, organic solvents.

These 1ml Inserts have a wide conical point for fast sample filling to minimize bubble formation and has outer dimensions of 7.72mm x 44.85mm.

Set includes [Pre-Racked Inserts](#), [Base](#) & [Pre-Slit Sealing Mat](#).

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$155.66**

Printed On 2024/05/20