


Catalog Number

95296-1mL-CS

Inserts, in Racks. Includes 96, Clear, 1000ul Glass Low Volume Precision Inserts with wide conical points. For use in the Micro-Sample Management System. U-2D Brand. 12/CS.



Additional Info

Inserts are 1st Hydrolytic, Class A, TYPE 33 Borosilicate Glass made for use with the U-2D Micro-Sample Management System. These 1ml Inserts have a wide conical point for fast sample filling to minimize bubble formation and has outer dimensions of 7.72mm x 44.85mm. Preloaded in Racks.

Click [HERE](#) for smaller pack quantity. *Bases, Sealing Film or Mats are sold separately.*

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$1,456.90

Printed On 2024/05/20