

Catalog Number

95296-LG-30P-EA

Mat, Sealing, Large, Silicone Rubber, Pre-Slit, Plug Style for 96, 500 and 1000ul Inserts. U-2D Brand. 1 EA.



Additional Info

Mats are made of ultra pure Silicone Rubber, designed and molded to seal up to 96 each 500ul and 1ml Inserts in the U-2D Micro-Sample Management System. These mats are Pre-Slit and have alphanumeric markings on the top side of the mats, along the edges for identification of the Inserts. Use with Aqueous Solutions.

Click [HERE](#) for case quantity.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$26.09

Printed On 2024/05/19