

Catalog Number

9532S-E-0AV

Vials, Screw Top, Glass. Amber / Clear Insert, 300ul Fused Insert, 9-425mm threads, 12x32mm outer dimensions. For use as an autosampler vial. BASIK Brand. 100/PK.



Additional Info

Vials and inserts are first hydrolytic, Class A, Borosilicate Glass, with a wide bottom, 300ul, inserts fused to the bottom of the vial. The wide bottom helps prevent bubble formation when filling quickly or with automation.

Our most economical Fused-Insert Vials. For budget minded laboratories.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$141.09

Printed On 2024/05/20