

## Catalog Number

9533S-E-0AV

**Vials, Screw Top, Glass. Amber / Clear Insert, 300ul Fused Insert™, 10-425mm threads, 12x32mm outer dimensions. For use as an autosampler vial. BASIK Brand. 100/PK.**



### Additional Info

Wide mouth vials and inserts are first hydrolytic, Class A, Borosilicate Glass, 300ul total fill volume, inserts fused to the bottom of the vial. Amber vials are often used for light sensitive compounds and are considered low actinic.

Our most economical Fused-Insert Vials with 10mm screw threads. For budget minded laboratories.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$280.35**

Printed On 2024/05/20