


Catalog Number

9582S-10X-K

Septa, 8mm, White Silicone Rubber / Red PTFE that is 0.045" thick and 0.321" diameter for use in 8-425mm thread caps. Caps sold separately. MicroSolv Brand. 1000/CS.



Additional Info

Used for both HPLC or GC, 8mm caps.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$192.50

Printed On 2024/05/20