

## Catalog Number

9582S-CAP1

**Caps, Screw Top**, with an Open Hole. "Knurled", polypropylene, Black Caps. Septa are sold separately and must be inserted manually. For use with 8-425mm thread, 12x32mm autosampler vials. BASIK Brand.



### Additional Info

Click [HERE](#) to visit a page with Septa for this Cap.

#### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

#### Price

**\$76.19**

Printed On 2024/05/20