

## Catalog Number

9582S-WAV

**Vials, Screw Top, Glass. Amber, 2ml, with Write-On Patch and fill lines, 8-425mm threads, 12x32mm outer dimensions. For use as an autosampler vial or storage bottle. BASIK Brand. 100/PK.**



### Additional Info

First hydrolytic, Class A, Borosilicate Glass, amberized. Write-On Patch and fill lines are epoxy paint which will resist smearing during use. Features a narrow mouth which may be useful to reduce evaporation and requires an 8mm screw cap.

Amber vials are often used for light sensitive compounds and considered low actinic. The narrow opening and write-on patch make these ideal for stability testing and sample storage.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$39.10**

Printed On 2024/05/20