

Catalog Number

9582S-WCV

Vials, Screw Top, Glass. Clear, 2ml, with Write-On Patch and fill lines, 8-425mm threads, 12x32mm outer dimensions. For use as an autosampler vial or storage bottle. BASIK Brand. 100/PK.



Additional Info

First hydrolytic, Class A, Borosilicate Glass. Write-On Patch and fill lines are epoxy paint which will resist smearing during use. Features a narrow mouth which may be useful to reduce evaporation.

The narrow opening and write-on patch make these ideal for stability testing and sample storage.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$33.90

Printed On 2024/05/20