


Catalog Number

96010-10-1501

Septa, 10mm, White Silicone Rubber / Beige PTFE that is 1.5mm thick for use in 10-425mm thread caps. Caps sold separately. MicroSolv Brand. 1000/CS.



Additional Info

Use for HPLC or GC.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$645.84

Printed On 2024/05/19