

Catalog Number

97004-01

Vial & Cap Kit Includes 144 EA 4ml, Screw Top, Clear Storage Vials and 144 EA matching, Attached Screw Caps with a Rubber / PTFE Liner fitted in the Black, Knurled, Phenolic Caps. BASIK Brand. Easy Purchase



Additional Info

These 4ml Vials are made with 1st Hydrolytic, TYPE 33 Borosilicate Glass with 13-425mm threads that are supplied in egg crate style corrugated tray, boxes.

Closed, or solid, Phenolic caps are supplied and can tolerate higher temperatures than polypropylene and are harder.

Caps and Vials are supplied assembled.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$482.31

Printed On 2024/05/19