


Catalog Number

97020-SP

Caps, Screw Top, Solid, with Red-Orange Natural Rubber / Transparent TFE Liner, fitted in knurled polypropylene, White Caps. For use with ND20 thread, storage Vials or Bottles. MicroSolv Brand. 100/PK.



Additional Info

Natural Rubber / TFE Liner is 60° Shore A (hardness) and is 1.3mm thick. Cap can be used for use with Storage Vials or Bottles with 20-400mm screw threads.

Use with Storage Vials [97020-CV-20](#).

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$77.75

Printed On 2024/05/19