

## Catalog Number

97031-0CV

**Vials, EPA Style**, Glass. Clear, 30ml. 70x27.5mm outer dimensions. For use as an EPA test vial or storage bottle. ND24 threads. BASIK Brand. 100/PK.



### Additional Info

First hydrolytic, Class A, Type 33, Borosilicate Glass. Fill volume is approximately 31ml.

Use 24mm screw caps that have 24-414 threads, sold separately. *Lower cost for budget minded laboratories.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$79.31**

Printed On 2024/05/20