

Catalog Number

97040-05C-K

Caps, Screw Top, with White Silicone Rubber / White PTFE Septa, fitted in knurled polypropylene, Black Caps. For use with 24-414mm threaded EPA Style vials. BASIK Brand. 100/PK..



Additional Info

EPA quality septa can be used for EPA testing or for long term storage with any 30, 40 or 60 EPA Style vial.

Cap has a 12.5mm top hole.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$67.01

Printed On 2024/05/31