

## Catalog Number

98100-SP13-K

**Caps, Screw Top, Solid, with a Silicone Rubber / PTFE liner, in "knurled", polypropylene, Black Caps. For use with 13-425mm thread, 15x45mm autosampler vials. MicroSolv Brand 100/PK.**



### Additional Info

Turns 4ml, 15x45 autosampler vials into storage bottles. The Silicone Rubber / PTFE liner provides an inert surface with good sealing properties.

#### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

#### Price

**\$41.51**

Printed On 2024/05/20