

Catalog Number

PQA-FKT-01

Flow Cell Filling Syringe. Replacement. Use with Chemical Solutions HSQ or PQ Kits. MicroSolv brand. 1 EA.



Additional Info

This Qualification accessory is used to fill flow cells on HPLC instruments for Performance Qualification of the detectors.

The spring (included) reduces the chance of putting too much pressure on the flow cell and adds the solution slowly.

This Flow Cell Filling Syringe is supplied assembled and ready to use.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$187.90

Printed On 2024/05/14