

## Catalog Number

PQA-RTM-01

**Resolution Test Mix Solution.** For HPLC qualification. 30 ml. Chemical Solutions brand. 1 EA.



### Additional Info

Use to test your HPLC or UV-vis detector for precision and accuracy with the Chemical Solutions PQ or HSQ kits. This replacement solution is labelled as Rs Test Mix (Resolution Test) and is NIST traceable. Supplied in a 30 mL plastic bottle and leak proof cap.

Click [HERE](#) for specifications and concentration of this solution.

This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$285.03**

Printed On 2024/05/09