

Catalog Number

04026-C

Capillary, Bare Fused Silica, 25um ID by 10 meters. Simplus Brand. 1 EA.



Additional Info

This untreated, fused silica capillary is clad in polyimide and designed for use in Capillary Electrophoresis. Supplied on a plastic, non static spool for storage.

Manufactured with extreme precision, these capillaries have excellent concentricity throughout the entire length of the capillary. The outer diameter (OD) is 365um. Made in USA.

Click [HERE](#) for 1 meter length and [HERE](#) for 25 meter length.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$682.59

Printed On 2024/05/02