


Catalog Number

04050-C

Capillary, Bare Fused Silica, 50um ID by 1 meter. Simplus Brand. 1 EA.



Additional Info

This untreated, fused silica capillary is clad in polyimide and designed for use in Capillary Electrophoresis.

Manufactured with extreme precision, these capillaries have excellent concentricity throughout the entire length of the capillary. The outer diameter (OD) is 365um. Made in USA.

Click [HERE](#) for 10 meter length and [HERE](#) for 25 meter length.

Price
\$164.28

Additional Info

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$164.28