

Catalog Number

04650-LF

Capillary, Controlled, Low Flow, 50um ID by 1 meter. MICROSOLV Brand. 1 EA.



Additional Info

This fused silica capillary is clad with a polyimide outer coating. The inner capillary wall is covalently bonded with a polymer with Sulfonic Acid Groups to create a lower than normal negative charge environment that minimally conducts current and creates a much slower than normal Electro-Osmotic Flow (EOF). Ionized analytes will move to the anode or cathode respectively, in longer migration times than on untreated Bare Fused Silica Capillaries..

This Capillary is meant to be cut to your own desired length for your own Capillary Electrophoresis Instrument and is supplied without finished tips.

Price
\$247.18

Additional Info

The Outer Diameter (OD) of this capillary is 365um.

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$247.18