


## Catalog Number

07220-T1

**Tubing, Replacement,** for Hand Held Low Volume Pump, 365um ID.  
Polyethylene, 1 meter long. MicroSolv Brand. 1 EA.



### Additional Info

This plastic, low pressure Tubing is a Replacement accessory for the MicroSolv Hand Held, Low Volume Pump Syringe for 365um OD CE or CEC Capillaries.

Slips over the Needle and Capillary with a negative tolerance.

Click [HERE](#) for Hand Held Pump ordering information and [HERE](#) for Replacement Syringe.

This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$142.59**

Printed On 2024/05/08