

Catalog Number

40000-05P-2

HPLC Column, Silica-C, 4um, 2.1mm ID x 50mm Length, 100A. Cogent TYPE-C Brand. 1 EA.



Additional Info

The Silica-C™ HPLC Columns are made with [Cogent TYPE-C™ Silica](#), and provide unique Selectivity due to its surface chemistry (especially in Normal Phase HPLC) compared to standard silica, This Silica Hydride Phase is not modified and therefore has no end-capping or organic ligands.

Chosen for analyzing compounds that are not soluble in Water including but not limited to: lipids, peptides, proteins, oligonucleotides and other compounds that require non polar solvents.

- Faster Equilibration than other Silica Columns
- No need to control Moisture in Mobile Phase Solvents
- Robust Column Lifetimes

This short and narrow-bore column provides short run times but requires a pump that can precisely deliver, and a flow cell that made for, reduced flow rates with lower volumes. This

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$839.68

Printed On 2024/05/03