

# Catalog Number

40000-25P-3

**HPLC Column, Silica-C, 4um, 3.0mm ID x 150mm Length, 100A.**  
Cogent TYPE-C Brand. 1 EA. CUSTOM ORDER.



## Additional Info

The Silica-C™ HPLC Columns are made with [Cogent TYPE-C™ Silica](#), and provide unique Selectivity due to its surface chemistry (especially in Normal Phase HPLC) compared to standard silica, This Silica Hydride Phase is not modified and therefore has no end-capping or organic ligands.

Chosen for analyzing compounds that are not soluble in Water including but not limited to: lipids, peptides, proteins, oligonucleotides and other compounds that require non polar solvents.

**Price****\$1,167.17**

## Additional Info

- Faster Equilibration than other Silica Columns
- No need to control Moisture in Mobile Phase Solvents
- Robust Column Lifetimes

*This size column provides similar run times as a 4.6mm ID column, but uses much less mobile phase, with lower flow rates. No special equipment is needed for this column to save solvent use.*

This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$1,167.17**