

## Catalog Number

## 40000-HG6

**Guard Column Cartridge, Silica-C, Replacement Cartridge, 4.0mm ID x 10mm Length, 4um, 100A. Hichrom style, in an individual black case. Cogent TYPE-C Brand. 1 EA.**



### Additional Info

This Guard Column is designed to protect your analytical column from particulate and contaminants that will adsorb to your column. Packed in a Stainless Steel Tube, it has SS Frits. Labelled, and shipped in an individual, black, Delrin case to protect and keep it from drying out. Use with [Hichrom Style Guard Column Holder](#).

Sometimes, the individual cartridge is preferred over buying a pack of 5 of the same phases for when different "Guard Column" phases are preferred for pairing them with different analytical phases for unique separations.

*Use this guard column cartridge with a 3.0, 4.0 and 4.6mm ID analytical column.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**

**\$273.86**

Printed On 2024/08/06

## Additional Info

**MicroSolv Technology Corporation**  
9158 Industrial Blvd. NE, Leland, NC 28451  
tel. (732) 380-8900, fax (910) 769-9435  
Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)  
Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**  
**\$273.86**

Printed On 2024/08/06