

Catalog Number

40000-PG4

Guard Column Starter Set, Silica-C, complete with 5 each 4.0 mm ID x 10 mm long cartridges packed with 4 μ m, 100 A phase and 1 each guard column holder. Cogent TYPE-C Brand. 1 EA.

**Additional Info**

Everything you need to get started using the Cogent TYPE-C™ guard columns with your methods. This kit includes one Stainless Steel [Guard Cartridge Column Holder](#) that can be used with any Cogent TYPE-C™ Guard Cartridge.

Five cartridges, packed with Cogent Silica-C stationary phase are included which will normally be sufficient protection for the lifetime of an analytical column, based on average use. Each cartridge in this kit is shipped in individual, closeable Delrin cases to keep them from drying out or being damaged.

Price**\$1,424.39**

Additional Info

Use this Starter Kit with 3.0, 4.0, and 4.6 mm ID analytical columns.

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$1,424.39