

Catalog Number

40008-PG2-3

Guard Column cartridges, Bidentate C8, replacement, 2.0 mm ID x 10 mm long cartridges packed with 5 μ m, 300 A phase. Cogent TYPE-C Brand. 5/PK.



Additional Info

These Guard Columns are designed to protect your analytical columns from particulate and contaminants that will adsorb to your column. These columns are packed in Stainless Steel tubes with SS Frits. Labelled, and shipped in individual, black, delrin cases to protect and keep them from drying out. Use with [guard column holder](#) sold separately.

Use these guard columns with 2.1 mm ID analytical columns.

MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451
tel. (732) 380-8900, fax (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com

Price
\$870.21

Printed On 2024/08/06

Additional Info

MicroSolv Technology Corporation
9158 Industrial Blvd. NE, Leland, NC 28451
tel. (732) 380-8900, fax (910) 769-9435
Email: customers@mtc-usa.com
Website: www.mtc-usa.com

Price
\$870.21

Printed On 2024/08/06