

## Catalog Number

## 40008-PG5

**Guard Column cartridges, Bidentate C8**, replacement, 4.0 mm ID x 10 mm long cartridges packed with 4  $\mu$ m, 100 A phase. Cogent TYPE-C Brand. 5/PK.



### Additional Info

These Guard Columns are designed to protect your analytical columns from particulate and contaminants that will adsorb to your column. These columns are packed in Stainless Steel tubes with SS Frits. Labelled, and shipped in individual, black, delrin cases to protect and keep them from drying out. Use with [guard column holder](#) sold separately.

*Use this guard column cartridge with a 3.0, 4.0, and 4.6 mm ID analytical columns.*

**MicroSolv Technology Corporation**  
9158 Industrial Blvd. NE, Leland, NC 28451  
tel. (732) 380-8900, fax (910) 769-9435  
Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)  
Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**  
**\$878.14**

Printed On 2024/08/06

## Additional Info

**MicroSolv Technology Corporation**  
9158 Industrial Blvd. NE, Leland, NC 28451  
tel. (732) 380-8900, fax (910) 769-9435  
Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)  
Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**  
**\$878.14**

Printed On 2024/08/06