

## Catalog Number

40008-PG6-3

**Guard Column cartridge, Bidentate C8**, replacement, 4.0 mm ID x 10 mm long cartridge packed with 5 um, 300 A phase. Cogent TYPE-C Brand. 1 EA.



### Additional Info

This Guard Column is designed to protect your analytical column from particulate and contaminants that will adsorb to your column. Packed in a Stainless Steel Tube, it has SS Frits. Labelled, and shipped in an individual, black, delrin case to protect and keep it from drying out. Use with [guard column holder](#) sold separately.

Sometimes, the individual cartridge is preferred over buying a pack of 5 when different "Guard Column" phases are preferred and pairing them with different analytical phases for unique separations.

*Use this guard column cartridge with a 3.0, 4.0, and 4.6 mm ID analytical columns.*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

**Price**

**\$236.58**

Printed On 2024/02/25