

## Catalog Number

40008-SP150

**Semi Prep HPLC Column, Bidentate C8, 4um, 100A. 10mm ID x 150mm Length. Cogent TYPE-C Brand. 1 EA.**



### Additional Info

This Cogent Bidentate C8 "Semi-Preparative", HPLC Column is made with 4um, 100A, [Cogent TYPE-Câ„¢ Silica](#), and this silica hydride based stationary phase provides excellent Selectivity, Faster Equilibration and Very Robust Column Lifetimes.

Chosen for Purification or "Prep" of Hydrophobic Molecules, Metabolites, and other Analytes with amine groups when a large volume of sample is desired to be collected when C18 is too retentive.

**Notes** - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$3,191.85**

Printed On 2024/07/27