

## Catalog Number

40018-P21-100

**Prep HPLC Column, Bidentate C18, 4um, 100A. 21.2mm ID x 100mm Length. Cogent TYPE-C Brand. 1 EA.**



### Additional Info

This Cogent Bidentate C18 Preparative, HPLC Column is made with 4um, 100A, [Cogent TYPE-C™ Silica](#), and this silica hydride based stationary phase provides excellent Selectivity, Faster Equilibration and Very Robust Column Lifetimes.

Chosen for Purification or "Prep" of Hydrophobic Molecules, Metabolites, and other Analytes with amine groups when a large volume of sample is desired to be collected.

**Notes** - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an

## Additional Info

order has been placed. Click [HERE](#) for more information on what is a Custom Order.

### **MicroSolv**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)