

Catalog Number

40018-SP50

Semi Prep HPLC Column, Bidentate C18, 4um, 100A. 10mm ID x 50mm Length. Cogent TYPE-C Brand. 1 EA.



Additional Info

This Cogent Bidentate C18 "Semi-Preparative", HPLC Column is made with 4um, 100A, [Cogent TYPE-Câ„¢ Silica](#), and this silica hydride based stationary phase provides excellent Selectivity, Faster Equilibration and Very Robust Column Lifetimes.

Chosen for Purification or "Prep" of Hydrophobic Molecules, Metabolites, and other Analytes with amine groups when a large volume of sample is desired to be collected.

Notes - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$2,913.29

Printed On 2024/07/27