

Catalog Number

40031-HG5

**Guard Column Cartridge, UDA (wcx), Replacement Cartridge, 2.0mm ID x 10mm Length, 4um, 100A. In an individual black case. Cogent TYPE-C Brand. 1 EA.**



### Additional Info

This Guard Column is designed to protect your analytical column from particulate and contaminants that will adsorb to your column. Packed in a Stainless Steel Tube, it has SS Frits. Labeled, and shipped in an individual, black, Delrin case to protect and keep it from drying out. Use with [Hichrom Style Guard Column Holder](#).

Sometimes, the individual cartridge is preferred over buying a pack of 5 when different "Guard Column" phases are preferred and pairing them with different analytical phases for unique separations.

Price

\$280.02

## Additional Info

*Use this guard column cartridge with a 2.1mm ID analytical column.*

### **MicroSolv**

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### **Price**

**\$280.02**