

Catalog Number

40036-15P

HPLC Column, Amide, 4um, 4.6mm ID x 150mm Length, 100A. Cogent TYPE-C Brand. 1 EA.



Additional Info

This Cogent Amide HPLC Column is made with [Cogent TYPE-C™ Silica](#), and this silica hydride based stationary phase provides excellent Selectivity, Faster Equilibration and Robust Column Lifetimes.

Chosen for analysis of Sugars and other Analytes with amine groups,

This column is one of the more popular column sizes and often chosen for use with standard LC systems that have ELSD or Corona CAD detectors or are working with UV absorbing nitrogen containing compounds.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$1,205.89

Printed On 2024/05/07