

Catalog Number

40036-P21-50

Prep HPLC Column, Amide, 4um, 100A. 21.2mm ID x 50mm Length.
Cogent TYPE-C Brand. 1 EA.



Additional Info

This Cogent Amide preparative, HPLC column is made with 4um, 100A, [Cogent TYPE-C™ Silica](#), and this silica hydride based stationary phase provides excellent selectivity, faster equilibration and robust column lifetimes.

Chosen for purification or "Prep" of sugars and other analytes with amine groups when a large volume of sample is desired to be collected.

Notes - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$6,247.84

Printed On 2024/05/03