

Catalog Number

40036-SP50

Semi Prep HPLC Column, Amide, 4um, 100A. 10mm ID x 50mm Length. Cogent TYPE-C Brand. 1 EA.



Additional Info

This Cogent Amide "Semi" Preparative, HPLC Column is made with 4um, 100A, Cogent [TYPE-C™ Silica](#), and this silica hydride based stationary phase provides excellent Selectivity, Faster Equilibration and Robust Column Lifetimes.

Chosen for Purification or "Prep" of Sugars and other Analytes with amine groups when a sample is desired to be collected.

Notes - A specialized HPLC System is required for columns this size. This is a CUSTOM ORDER manufactured item which cannot be returned or canceled once an order has been placed. Click [HERE](#) for more information on what is a Custom Order.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$2,913.29

Printed On 2024/05/05