

## Catalog Number

40100-15P-2

**Method Development Kit 2154, TYPE-C.** Includes 5 Cogent HPLC Columns, 2.1mm x 150mm. Cogent TYPE-C Brand. 1 PK.



### Additional Info

This Kit, 2154 is an economical way to purchase 5 Cogent HPLC Columns that are often used with LCMS for Method Development and Scouting. Each Column has a 2.1mm Internal Diameter and is 150mm Long, packed with 4um particles, 100A pore size and are manufactured with [Cogent TYPE-C™ Silica Hydride](#).

This kit includes the columns listed below:

- **Cogent Diamond Hydride™**, Click [HERE](#) for More Information.
- **Cogent Bidentate C18™**, Click [HERE](#) for More Information.
- **Cogent Bidentate C8™**, Click [HERE](#) for More Information.
- **Cogent Silica-C™**, Click [HERE](#) for More Information.
- **Cogent UDC-Cholesterol™**, Click [HERE](#) for More Information.

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$5,518.36**

Printed On 2024/05/04