

Catalog Number

40218-15P-2

HPLC Column, Bidentate C18 2.o, 2.2um, 2.1mm ID x 150mm Length, 120A. Cogent TYPE-C Brand. 1 EA.



Additional Info

This Cogent Bidentate C18 2.o HPLC Column is made with [Cogent TYPE-C™ Silica](#), and this smaller particle, silica hydride based stationary phase provides Increased Efficiency, Excellent Selectivity and Faster Equilibration. Considered "Near UHPLC", this column can be used on any HPLC Instrument and still deliver improved efficiency.

Chosen for analysis of hydrophobic compounds that are difficult to retain in Aqueous Mobile Phases, including very non-polar compounds.

This column can be used for LCMS when increased resolution is desired but an LC

Price**\$1,503.14**

Additional Info

instrument designed for columns of this ID is required.

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$1,503.14