

## Catalog Number

40231-HG3

**Guard Column Cartridge, UDA (wcx) 2.o**, Replacement Cartridge, 2.0mm ID x 10mm Length, 2.2um, 120A. Hichrom style, in an individual black case. Cogent TYPE-C Brand. 1 EA.



### Additional Info

This Guard Column is designed to protect your analytical column from particulate and contaminants that will adsorb to your UDA (wcx) 2.o Column. Packed in a Stainless Steel Tube, it has SS Frits. Labelled, and shipped in an individual, black, Delrin case to protect and keep it from drying out. Use with [Hichrom Style Guard Column Holder](#).

Sometimes, the individual cartridge is preferred over buying a pack of 5 when different "Guard Column" phases are preferred and pairing them with different analytical phases for unique separations.

*Use this guard column cartridge with a 2.1mm ID analytical column packed with Undecanoic Acid 2.o (UDA).*

### MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$275.19**

Printed On 2024/05/02