

Catalog Number

40231-PG2

Guard Column cartridges, UDA 2.0 (wcx), replacement, 2.0 mm ID x 10 mm long cartridges packed with 2.2 μ m, 120 A phase. Cogent TYPE-C Brand. 5/PK.



Additional Info

These Guard Columns are designed to protect your analytical columns from particulate and contaminants that will adsorb to your column. These columns are packed in Stainless Steel tubes with SS Frits. Labelled, and shipped in individual, black, Delrin cases to protect and keep them from drying out. Use with [guard column holder](#) sold separately.

Use these guard columns with 2.1 mm ID analytical columns.

Price

\$1,043.96

Additional Info

MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$1,043.96