

Catalog Number

43000-15P

Metabolomics Development Kit, 4.6x150. 1 Cogent Diamond Hydride™ & 1 Cogent Bidentate C18™. Cogent TYPE-C Brand. 1 PK.



Additional Info

This is an economical way to purchase two HPLC Columns that are often used in Metabolomics and many other fields.

Each Column has a 4.6mm Internal Diameter and is 150mm Long, packed with 4um particles and 100A pore size and are manufactured with [Cogent TYPE-C™ Silica Hydride](#).

- Click [HERE](#) for More Information about Cogent Diamond Hydride™
- Click [HERE](#) for More Information about Cogent Bidentate C18™.

MicroSolv Technology Corporation

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: customers@mtc-usa.com

Website: www.mtc-usa.com

Price

\$1,992.53

Printed On 2024/05/03