

# Catalog Number

43200-05P-2

**Metabolomics Development Kit, 2.o, 2x50.** 1 Cogent Diamond Hydride 2.o™ & 1 Cogent Bidentate C18 2.o™ 2.1mm 50mm. Cogent TYPE-C Brand. 1 PK.



## Additional Info

This is an economical way to purchase two HPLC Columns that are often used with LCMS in Metabolomics and many other fields.

Each Column has a 2.1mm Internal Diameter and is 50mm Long, packed with 2.2um particles and 120A pore size and are manufactured with [Cogent TYPE-C™ Silica Hydride](#).

- Click [HERE](#) for More Information about Cogent Diamond Hydride™
- Click [HERE](#) for More Information about Cogent Bidentate C18™.

**Price****\$2,270.54**

## Additional Info

### MicroSolv

9158 Industrial Blvd. NE, Leland, NC 28451

tel. (732) 380-8900, fax (910) 769-9435

Email: [customers@mtc-usa.com](mailto:customers@mtc-usa.com)

Website: [www.mtc-usa.com](http://www.mtc-usa.com)

### Price

**\$2,270.54**